

# DARCO<sup>®</sup>

## Headed Compression Screw



# DARCO® Headed Compression Screw

Versatile cannulated screw designed specifically for foot and ankle applications

- » Full range of screws diameters:  
Available in 2.7, 3.5, 4.0, 4.5, 6.5, 7.5mm
- » Titanium-alloy construction
- » Available in various thread lengths for ultimate surgical flexibility
- » Short, long and fully threaded options available



## Alternative Products:

CHARLOTTE® Multi-Use Compression Screws  
DARCO® Headless Compression Screws

## Complementary Products:


ORTHOLOC® MDCO Plates  
TCP Total Compression Plates  
CHARLOTTE® CLAW®II Plates  
CHARLOTTE® Snap-Off Screws  
CHARLOTTE® Quick Staple



**ORTHOLOC® 3Di**  
MDCO Plate



**TCP**  
Total Compression Plate



**CHARLOTTE®**  
Snap-Off Screw



**CHARLOTTE® CLAW®II**  
Polyaxial Compression  
Plating System



1023 Cherry Road  
Memphis, TN 38117  
800 238 7117  
901 867 9971  
www.wright.com

62 Quai Charles de Gaulle  
69006 Lyon  
France  
+33 (0)4 72 84 10 30  
www.tornier.com

Unit 1, Campus Five  
Letchworth Garden City  
Hertfordshire SG6 2JF  
United Kingdom  
+011 44 (0)845 833 4435